




Hernia Repair and Soft Tissue Augmentation with Strength,
Softness & Sustainability

FILAPROP™ MESH Soft

Non Absorbable Soft Polypropylene Mesh



Indications

- ◆ Open & Laparoscopic Inguinal Hernia
- ◆ Open Ventral Hernia Repair
- ◆ Fascial Deficiencies
- ◆ Chest Wall Reconstruction




Key Features

- 60% lighter weight than conventional polypropylene mesh.
- Affords excellent strength, durability & fast tissue integration.
- Soft & pliable mesh with unique knitting which facilitates natural process of wound healing.

Technical Specifications

Pore Size	0.8 mm x 1.0 mm
Weight	45 gsm
Thickness	0.34 mm
Burst Strength	81.29 N/cm
Sterilization Method	Ethylene Oxide

Product Specifications



Filaprop™ Mesh Soft Non-Absorbable Surgical Soft Mesh

Sr. No.	Product Code	Brand Name	Description
1	PPMS611	FILAPROP MESH SOFT	POLYPROPYLENE MESH SOFT (06 CM X 11 CM)
2	PPMS715	FILAPROP MESH SOFT	POLYPROPYLENE MESH SOFT (7.6 CM X 15 CM)
3	PPMS1015	FILAPROP MESH SOFT	POLYPROPYLENE MESH SOFT (10 CM X 15 CM)
4	PPMS1215	FILAPROP MESH SOFT	POLYPROPYLENE MESH SOFT (12 CM X 15 CM)
5	PPMS1218	FILAPROP MESH SOFT	POLYPROPYLENE MESH SOFT (12 CM X 18 CM)
6	PPMS1515	FILAPROP MESH SOFT	POLYPROPYLENE MESH SOFT (15 CM X 15 CM)
7	PPMS1520	FILAPROP MESH SOFT	POLYPROPYLENE MESH SOFT (15 CM X 20 CM)
8	PPMS1530	FILAPROP MESH SOFT	POLYPROPYLENE MESH SOFT (15 CM X 30 CM)
9	PPMS3030	FILAPROP MESH SOFT	POLYPROPYLENE MESH SOFT (30 CM X 30 CM)