

MIRUS™ Skin Stapler

Outstanding performance across wider skin wounds

Skin staplers are typically used at the end of a surgical procedure for routine wound closure. The staples are then removed with staple extractor after the skin has had an opportunity to heal.

The MIRUS™ Skin Stapler's handle design is characterized by good texture for an user friendly experience.

The angled head provides a clear view to ensure precise staple placement and allows the pins to staple the tissue easily.

Staple material : 316L Stainless Steel

Sterilization : Ethylene Oxide

Packaging : 8 pcs/Box


Code : MSS35

Size : 35

Staples : Regular



Staple Dimension:



Larger diameter wire, greater width, and extra staple height leads to a secure wound closure

Stapler	Skin Stapler
Procedure Application	Wound closure